BIOHIT

Innovating for Health

Biohit Oyj

ColonView Hb/Hp Test

Instructions for use

Käyttöohje

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Istruzioni d'uso

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Biohit Oyj
Laippatie 1
FIN-00880 Helsinki
Finland
Tel: +358-9-773 861
Fax: +358-9-773 86200
E-mail: info@biohit.com
www.biohit.com

REF 602 200
602 210
602 230

IVD

CE

EXPLANATION OF SYMBOLS USED
KÄYTETTYJEN SYMBOLIEN SELITYKSET
FÖRKLARING TILL DE ANVÄNDA SYMBOLERNA
ERKLÄRUNG DER SYMBOLE AUF ETIKETTEN

	IVD	REF	LOT
ENGLISH	For in vitro diagnostic use	Catalogue number	Batch code
SUOMI	In vitro diagnostiseen käyttöön	Kataloginnumero	Eränumero
SVENSKA	För diagnostisk användning in vitro	Katalog numret	Batch kod
DEUTSCH	Zur Verwendung für die In-vitro-Diagnose	Artikelnummer	Chargennummer
FRANCAIS	Pour le diagnostique in vitro	Numéro de Catalogue	Code du lot
ITALIANO	Per l'uso diagnostico in vitro	Numero di catalogo	Codice lotto
ESPAÑOL	Para uso diagnóstico in vitro	Número de catálogo	Código de lote

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CONTENTS

1. INTENDED USE.....	6
2. CLINICAL BACKGROUND.....	6
3. PRINCIPLE OF THE COLONVIEW HB/HP TEST	7
4. KIT CONTENTS AND REAGENT PREPARATION	7
4.1 Test Cassettes	7
4.2 ColonView Sample Collection Kits	8
4.3 Instructions for Use (x 1)	8
5. MATERIALS REQUIRED BUT NOT PROVIDED	8
6. STORAGE AND STABILITY	8
7. WARNINGS AND PRECAUTIONS	8
8. SPECIMEN COLLECTION AND HANDLING.....	9
9. TEST PROCEDURE	10
9.1 Procedure	10
9.2 Quality Control/ Internal Procedural Control.....	11
10. RESULTS	11
11. LIMITATIONS OF THE PROCEDURE.....	12
12. PERFORMANCE CHARACTERISTICS.....	13
13. DATE OF ISSUE	14
14. WARRANTY	14

1. INTENDED USE

ColonView Hb/Hp Fecal Occult Blood Test is a visual, immuno-chromatography quick test for the qualitative detection of human hemoglobin/haptoglobin complex (Hb/Hp) in stool samples. For professional, *in vitro* diagnostic use.

2. CLINICAL BACKGROUND

ColonView Hb/Hp Fecal Occult Blood Test (FOB) test was designed to aid diagnosis of lower GI pathologies, such as colorectal cancers and large adenomas that bleed. Colorectal cancer is one of the most commonly diagnosed cancers and a leading cause of cancer death in the United States (1-2). Screening for colorectal cancer probably increases the cancer detection rate at an early stage, therefore reducing the mortality (3-5).

Earlier commercially available FOB tests utilized the guaiac test, which requires special dietary restriction to minimize false positive and false negative results. ColonView Hb/Hp tests are specially designed to detect human Hb/Hp in stool samples using immunochemical methods, which improved specificity for the detection of lower GI disorders, including colorectal cancers and adenomas (6-10).

A molecule of Hb consists of 2 pairs of peptide (alpha and beta globins) chains and four heme groups, each with one atom of iron. Free Hb may separate into alpha-beta molecules, which are bound to a protein called haptoglobin. Hb/Hp complex plays an important role in the retrieval of hemoglobin from lysed erythrocytes and is relatively stable against acid and proteolytic degradation. Therefore the Hb/Hp complex can be detected even after longer passages through the bowel. This means that blood caused by larger intestinal polyps and higher carcinomas can also be detected. Detection of the Hb/Hp complex displays a significant in-

crease in sensitivity with regard to the recognition of colorectal adenomas, and with regard to also carcinomas when combined with Hb detection (8-9).

3. PRINCIPLE OF THE COLONVIEW HB/HP TEST

The immunological ColonView Hb/Hp test is based on an immunochromatographic method, in which Hb/Hp is specifically recognized through specific antibody reactions.

In the case of a positive result, Hb/Hp molecules from the stool sample attach to the test band (T) and are visible by means of a pink/red coloration. In the case of a negative result, there are no Hb/Hp molecules that can attach to the test band (T), and therefore there is no coloration of the test band (T). If the control strip (C) turns red/pink in color, this shows that the sample has been taken correctly and has migrated correctly. The test is therefore valid.

4. KIT CONTENTS AND REAGENT PREPARATION

Each ColonView Hb/Hp test pack contains reagents for 10 (Cat. No. 602200, 602210) or 30 (Cat.No. 602230) tests. Each test pack contains:

4.1 Test Cassettes

Contents: Cat.No. 60200: 10 individually wrapped Test Cassettes without Sample Collection Kits.

Cat.No.602210: 10 individually wrapped Test Cassettes.

Cat.No. 602230: 30 individually wrapped Test Cassettes.

Preparation: Ready to use.

4.2 Sample Collection Kits

Contents: Cat.No. 602210: 10 x ColonView Sample Collection Kit (1). Each kit contains one Sample Collection Tube, one Stool Collection Paper and one Instruction for Sample Collection.

Cat.No. 602230: 10 x ColonView Sample Collection Kit (3). Each kit contains three Sample Collection Tubes, three Stool Collection Papers and one Instruction for Sample Collection.

Sample Collection Tubes contain 2 ml of sample buffer with BSA, 0.09% sodium azide and preservative.

Preparation: Ready to use

4.3 Instructions for Use (x 1)

5. MATERIALS REQUIRED BUT NOT PROVIDED

A piece of tissue paper to prevent solution from splashing when opening the Sample Collection Tube.

6. STORAGE AND STABILITY

The Test Cassettes and ColonView Sample Collection Kits can be stored at +2° - +30°C until the expiration date.

7. WARNINGS AND PRECAUTIONS

- Valid only for *in vitro* diagnostic use and only for professional use.

- Handle stool specimens as potential biohazardous material.
- It is recommended that the Stool Collection Paper enclosed is used, in order to avoid mixing with blood from urine samples or with interfering components from the toilet water.
- The saline solution in the Sample Collection Tube contains a small amount of sodium azide, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of buffered saline or extracted samples, always flush with copious quantities of water to prevent azide build up.
- Do not use the Sample Collection Tube or Test Cassette after the expiration date.
- Only open the Test Cassette packaging if the test is actually going to be performed.
- Do not spill solution onto the Test Cassette apart from placing sample into the sample well.
- Patients should follow the instructions for obtaining a sample exactly and should not take any sample during menstruation or in case of bleeding hemorrhoids.
- Bring all reagents to room temperature (20-30°C) at least 10 minutes before use.
- Evaluate the test result after 15 minutes and not beyond 25 minutes.

8. SPECIMEN COLLECTION AND HANDLING

Collection of 1 stool sample: follow instructions below. Collection of 3 stool samples: Collect the samples from 3 different stools over no more than 3 consecutive days.

1. Collect a specimen of stool in a clean dry container or receptacle, for example the Stool Collection Paper provided

with the Sample Collection Kit. Use the Stool Collection Paper as follows:

- Peel off the liner covering the adhesive tape on each end of the collection paper.
 - Lift toilet seat. Unfold collection paper and place on the rim of toilet bowl. Secure the adhesive tabs on the collection paper to the sides of the toilet rim.
 - Make bowel movement onto collection paper.
2. Unscrew (blue end) and remove the Sample Collection Tube applicator stick. Be careful not to spill or spatter solution from container.
 3. Collect random samples by inserting the applicator stick into the stool specimen. Take samples from various surfaces of the stool specimen.
 4. Re-insert the applicator stick into the tube and screw the cap tightly. Be careful not to break the tip of the Sample Collection Tube (white end).
 - If you are using the Stool Collection Paper, release adhesive tabs and flush the collection paper with stool.
 5. Write your name and the date of sample collection in the space provided, on the Sample Collection Tube.
 6. Return the specimen for testing promptly by mail or in person. The specimen(s) can be stored in the refrigerator (2-8 °C) for no more than 11 days, or at room temperature (max. 25 °C) for no more than 5 days.

9. TEST PROCEDURE

9.1 Procedure

1. The Test Cassette and the Sample Collection Tube containing the stool sample should be brought to room temperature

(+20°C to +30°C) at least 10 minutes before testing.

2. Take the required number of test cassettes from the foil packaging only immediately before performing the test. Mark the cassette with the name of the patient or with another form of identification.
3. Carefully shake the Sample Collection Tube to ensure that the stool sample mixes properly with the saline solution.
4. Unscrew the protective white cap. Take a paper towel and break the seal of the Sample Collection Tube with rotary motion. Hold the collection container upright and slowly add 3 drops of the solution into the round sample window (S) of the Test Cassette.
5. Read the results after 15 minutes. Strongly positive results may be evaluated even sooner. Evaluate the result a maximum of 25 minutes after adding the solution into the round sample window.

9.2 Quality Control/ Internal Procedural Control

The test contains a procedural control. A colored line that appears in the control region (C) shows that each test is performed correctly. It is easy for the background to become slightly yellowish in color during testing, depending on the color of the stool sample. This is acceptable, as long as evaluation of the test results is not adversely affected. The test is invalid if the result lines are not visible.

10. RESULTS

Test result is evaluated as “positive” if two lines appear.


Positive: 2 red/pink-colored lines appear; one in the control region

(C) and one in the test region (T).

If you are testing strongly positive samples, the intensity of the control line may be reduced. It is not recommended to compare the intensity of the lines.

Negative: 1 red/pink line appears in the control region (C).

Invalid: If no red/pink line appears in control region (C), this is a sign that the test is not functioning properly, or that the test materials are not correct. In this case, repeat the test with a new Test Cassette or contact the manufacturer for technical support.

11. LIMITATIONS OF THE PROCEDURE

A positive test occurs if human hemoglobin/haptoglobin complex is present in the sample.

Not all intestinal bleeding is the result of benign or malignant polyps. The information that you obtain by performing this test should be examined in conjunction with other clinical findings and test methods.

Because carcinomas and polyps can bleed to different extents and intermittently, it is advisable to test various stool samples from

eg. three different days.

Urine and excessive dilution of samples with water from the toilet bowl can lead to false test results. It is therefore recommended that the Stool Collection Paper be used.

Stool samples should not be taken during menstruation or 3 days before or after, in the case of bleeding caused by constipation, bleeding hemorrhoids, or in the case of medicine being administered rectally. This could lead to false positive results.

Blood-thinning medications such as, for example, acetylsalicylic acid or coumarin, as well as iron preparations, can lead to bleeding not caused by tumors.

12. PERFORMANCE CHARACTERISTICS

Analytical sensitivity

As low as 25 µg of hemoglobin/haptoglobin per liter of buffer solution produces a positive result. Prozone effect: The test also works reliably at extremely high Hb/Hp values (>500 mg/L).

Specificity

The Biohit ColonView Hb/Hp test is specific to human haemoglobin/haptoglobin complex and has no cross-reactivity with blood from cattle, pigs, rabbits or poultry. The Hb/Hp test does cross-react with blood from goat.

Note: In a study of the Shinshu University School of Medicine in Japan (11) the cost-value ratio of multiple measurements was examined. It shows the relative sensitivity increases with the number of consequent tests, and its relative specificity decreases slightly (see below).

Number of tests	sensitivity	specificity
1	58%	96%
2	89%	95%
3	100%	94%

13. DATE OF ISSUE

ColonView Hb/Hp test Instructions for use.
Version 02, 08.10.2008.

14. WARRANTY

The Manufacturer shall remedy all defects discovered in any Product (the “Defective Product”) that result from unsuitable materials or negligent workmanship and which prevent the mechanical functioning or intended use of the Products including, but not limited to, the functions specified in the Manufacturer’s specifications for the Products. ANY WARRANTY WILL, HOWEVER, BE DEEMED AS VOID IF FAULT IS FOUND TO HAVE BEEN CAUSED BY MALTREATMENT, MISUSE, ACCIDENTAL DAMAGE, INCORRECT STORAGE OR USE OF THE PRODUCTS FOR OPERATIONS OUTSIDE THEIR SPECIFIED LIMITATIONS OR OUTSIDE THEIR SPECIFICATIONS, CONTRARY TO THE INSTRUCTIONS GIVEN IN THE INSTRUCTION MANUAL.

The period of this warranty for the Distributor is defined in the instruction manual of the Products and will commence from the date the relevant Product is shipped by the Manufacturer. In case of interpretation disputes the English version applies. All Biohit Diagnostic kits have been manufactured according to our ISO 9001 / ISO 13485 quality management protocols and have passed all relevant Quality Assurance procedures related to these products.

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1. Lieberman D. Screening/Early Detection Model for Colorectal Cancer, Why Screen? *Cancer Supplement* 1994; 74(7): 2023-2027.
2. *Morbidity and Mortality Weekly Report*; (2007) 56(10): 227-228.
3. Dam, JV et. Al. Fecal Blood Screening for Colorectal Cancer. *Archive of Internal Medicine* 1995;155:2389-2402
4. Miller AB. An Epidemiological Perspective on Cancer Screening. *Clinical Biochemistry* 1995; 28(1): 41-48.
5. Ransohoff DF and Lang CA. Improving the Fecal Occult- Blood Test. *The New England Journal of Medicine* 1996; 334 (3): 189-190.
6. Frommer DJ et al. Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood. *British Medical Journal* 1988; 296:1092-1094
7. St. John DJB et al. Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia. *Gastroenterology* 1993; 104:1661-1668
8. Lüthgens K et al. Hemoglobin-Haptoglobin Complex: A Highly Sensitive Assay for the Detection of Fecal Occult Blood. *Clinical Laboratory* 1998; 44: 543-551.
9. Sieg A et al. Detection of colorectal neoplasms by the highly sensitive hemoglobin-haptoglobin complex in feces. *International Journal of Colorectal Disease* 1999; 14: 267-271.
10. van Rossum L et al. Random comparison of guaiac and immunochemical fecal occult blood tests for colorectal cancer in screening population. *Gastroenterology* 2008; 135: 82-90.
11. Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods : *Hepatogastroenterology*; 2000 Mar-Apr; 47(32):396-9.

**ORDERING INFORMATION/TILAUSTIEDOT/ BESTÄLLNINGSS
INFORMATION/ BESTELLINFORMATIONEN/ INFORMATIONS
RELATIVES AUX COMMANDES/ INFORMAZIONI RELATIVE
ALL'ORDINAZIONE/ INFORMACIÓN PARA LOS PEDIDOS**

ColonView Hb/Hp Test.

602200 (10 Test Cassettes only), 602210 (10), 602230 (30).

Headquarters
Biohit Oyj
Laippatie 1
00880 Helsinki, Finland
Tel: +358-9-773 861
Fax: +358-9-773 86290
E-mail: info@biohit.com
www.biohit.com

U.S.A.
Biohit. Inc.
Tel: +1-732-922-4900
Fax: +1-732-922-0557
E-mail:
customersvc.usa@biohit.com
biohit.com
(pipet@biohit.com)

CHINA
Finland Biohit Co., Ltd
Tel: +86-21-62485589
Fax: +86-21-62487786
E-mail: info.china@biohit.com

RUSSIA
Biohit OOO
Tel: +7-812-327 5327
Fax: +7-812-327 5323
E-mail: main@biohit.ru

FRANCE
Biohit S. A.
Tel: +33-1-3088 4130
Fax: +33-1-3088 4102
E-mail:
commercial.france@biohit.com

U. K.
Biohit Ltd.
Tel. +44-1803 315 900
Fax: +44-1803 315 530
E-mail: info@biohit.co.uk

GERMANY
Biohit Deutschland GmbH
Tel: +49-6003-82 820
Fax: +49-6003-828 222
E-mail: info@biohit.de

JAPAN
Biohit Japan Co., Ltd.
Tel: +81-3-5822 0021
Fax: +81-3-5822 0022
E-mail: sales@biohit.co.jp