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## BIBLIOGRAFIA

1. Ransohoff DF and Lang CA. Screening for colorectal cancer with the Fecal Occult Blood Test: a background paper. Ann Intern Med. 1997; 126: 811-822.
2. Towler BP, Irwig L, Glasziou P, Weller D, Kewenter J. Screening for colorectal cancer using the faecal occult blood test, Hemoccult. Cochrane Database Syst Rev. 2000;(2): CD001216.
3. Ransohoff DF and Lang CA. Suggested technique for Faecal Occult Blood testing and interpretation in colorectal cancer screening. Ann Intern Med. 1997; 126: 808-810.
4. Chien-Hua Chiang et al., "A Comparative Study of Three Fecal Occult Blood Test in Upper Gastrointestinal Bleeding". Kaohsiung J. Med. Sci. May 2006, 22; 223-228.
1. Virtanen et al., "Higher concentrations of serum transferrin receptor in children than in adults", Am. J. Clin. Nutr. 1999; 69: 256-60.
2. Angriman I. et al. Enzymes in feces: Useful markers of chronic inflammatory bowel disease. Clinica Chimica Acta 381 Feb 2007, p. 63-68
3. Quail, M.A. et al. Fecal Calprotectin Complements Routine Laboratory Investigations in Diagnosing Childhood Inflammatory Bowel Disease.

- Inflamm Bowel Dis, Vol 15 No 5; May 2009, p. 756-759
5. Gaya D.R., et al. Faecal calprotectin in the assessment of Crohn's disease activity. *Q J Med* 2005; Vol 98, May 2005, p. 435-441.
  6. Langhorst, M.D. et al. Noninvasive Markers in the Assessment of Intestinal Inflammation in Inflammatory Bowel Diseases: Performance of Faecal Lactoferrin, Calprotectin and PMN-Elastase, CRP, and Clinical Indices. *Am. J. Gastroenterol.* 2008; Vol 103, p. 162-169.
  8. Berntzen HB, Olmez U, Fagerhol MK, Munthe E. The leukocyte protein L1 in plasma and synovial fluid from patients with rheumatoid arthritis and osteoarthritis. *Scand J Rheumatol* 1991; 20(2): 74-82.
  9. Bjerke K, Halstensen TS, Jahnsen F, Pulford K, Brandtzaeg P. Distribution of macrophages and granulocytes expressing L1 protein (calprotectin) in human Peyer's patches compared with normal ileal lamina propria and mesenteric lymph nodes. *Gut* 1993; 34(10): 1357-63.
  10. Bjarnason I, Sherwood R. Faecal calprotectin: a significant step in the noninvasive assessment of intestinal inflammation. *J Pediatr Gastroenterol Nutr* 2001; 33: 11.
  11. Brandtzaeg P, Gabrielsen TO, Dale I, Muller F, Steinbakk M, Fagerhol MK. The leukocyte protein L1 (calprotectin): a putative nonspecific defence factor in epithelial surfaces. *Adv Exp Med Biol* 1995; 371A: 201-6.
  12. Bunn SK, Bisset WM, Main MJ, Golden BE. Faecal calprotectin as a measure of disease activity in childhood inflammatory bowel disease. *J Pediatr Gastroenterol Nutr* 2001; 32(2): 171-7
  13. Campeotto F, Baldassarre M, Butel MJ, Viallon V, Nganzali F, Soulaïnes P, Kalach N, Lapillonne A, Laforgia N, Moriette G, Dupont C, Kapel N. Faecal calprotectin: cutoff values for identifying intestinal distress in preterm infants. *J Pediatr Gastroenterol Nutr.* 2009 Apr;48(4):507-10.
  14. Campeotto F, Baldassarre M, Butel MJ, Viallon V, Nganzali F, Soulaïnes P, Kalach N, Lapillonne A, Laforgia N, Moriette G, Dupont C, Kapel N. Faecal calprotectin: cutoff values for identifying intestinal distress in preterm infants. *J Pediatr Gastroenterol Nutr.* 2009 Apr;48(4):507-10.
  15. Campeotto F, Kalach N, Lapillonne A, Butel MJ, Dupont C, Kapel N. Time course of faecal calprotectin in preterm newborns during the first month of life. *Acta Paediatr.* 2007 Oct;96(10):1531-3. Epub 2007 Aug 20. No abstract available.
  16. Campeotto F, Butel MJ, Kalach N, Derriex S, Aubert-Jacquin C, Barbot L, Francoal C, Dupont C, Kapel N. High faecal calprotectin concentrations in newborn infants. *Arch Dis Child Fetal Neonatal Ed.* 2004 Jul;89(4):F353-5.
  17. Canani RB, Terrin G, Rapacciuolo L, Miele E, Siani MC, Puzone C, Cosenza L, Staiano A, Troncone R. Faecal calprotectin as reliable non-invasive marker to assess the severity of mucosal inflammation in children with inflammatory bowel disease. *Dig Liver Dis.* 2008 Jul;40(7):547-53. Epub 2008 Mar 20.
  18. Canani RB, de Horatio LT, Terrin G, Romano MT, Miele E, Staiano A, Rapacciuolo L, Polito G, Bisesti V, Manguso F, Vallone G, Sodano A, Troncone R. Combined use of noninvasive tests is useful in the initial diagnostic approach to a child with suspected inflammatory bowel disease. *J Pediatr Gastroenterol Nutr.* 2006 Jan;42(1):9-15.
  19. Carroccio A, Iacono G, Cottone M, Di Prima L, Cartabellotta F, Cavataio F, Scalici C, Montalto G, Di Fede G, Rini G, Notarbartolo A, Averna MR. Diagnostic accuracy of faecal calprotectin assay in distinguishing organic causes of chronic diarrhea from irritable bowel syndrome: a prospective study in adults and children. *Clin Chem.* 2003 Jun;49(6 Pt 1):861-7.
  20. Carroccio A, Rocco P, Rabitti PG, Di Prima L, Forte GB, Cefalà AB, Pisello F, Geraci G, Uomo G. Plasma calprotectin levels in patients suffering from acute pancreatitis. *Dig Dis Sci.* 2006 Oct;51(10):1749-53.
  21. Costa F, Mumolo MG, Ceccarelli L, Bellini M, Romano MR, Sterpi C, Ricchiuti A, Marchi S, Bottai M. Calprotectin is a stronger predictive marker of relapse in ulcerative colitis than in Crohn's disease. *Gut.* 2005 Mar;54(3):364-8.
  22. Costa F, Mumolo MG, Bellini M, Romano MR, Ceccarelli L, Arpe P, Sterpi C, Marchi S, Maltinti G. Role of faecal calprotectin as non-invasive marker of intestinal inflammation. *Dig Liver Dis.* 2003 Sep;35(9):642-7.
  23. D'Inca R, Dal Pont E, Di Leo V, Benazzato L, Martinato M, Lamboglia F, Oliva L, Sturmiolo GC. Can calprotectin predict relapse risk in inflammatory bowel disease? *Am J Gastroenterol.* 2008 Aug;103(8):2007-14.
  24. D'Inca R, Dal Pont E, Di Leo V, Ferronato A, Fries W, Vettorato MG, Marines D, Sturmiolo G. Calprotectin and lactoferrin in the assessment of intestinal inflammation and organic disease. *Int J Colorectal Dis.* 2007 Apr;22(4):429-37. Epub 2006 Jul 13.
  25. Dale I, Fagerhol MK, Naesgaard I. Purification and partial characterization of a highly immunogenic human leukocyte protein, the L1 antigen. *Europ J Biochem* 1983; 134(1): 1-6.
  26. Fagerhol MK. Calprotectin, a faecal marker of organic gastrointestinal abnormality. *Lancet* 2000; 356(9244):1783-4.
  27. Fagerhol MK, Andersson KB, Naess-Andresen CF, Brandtzaeg P, Dale I. Calprotectin (The L1 Leukocyte Protein) In: VL Smith & JR Dedman (Eds): Stimulus Response Coupling, The Role of Intracellular Calcium-Binding Proteins. CRC Press, Boca Raton, Fla., USA, 1990, pp. 187-210.
  28. Gilbert JA, Ahlquist DA, Mahoney DW, Zinsmeister AR, Rubin J, Ellefson RD. Faecal marker variability in colorectal cancer: calprotectin versus hemoglobin. *Scand J Gastroenterol* 1996; 31(10):1001-5.
  29. John B, Kronborg O, Ton HI, Kristinsson J, Fuglerud P. A new faecal calprotectin test for colorectal neoplasia. Clinical results and comparison with previous method. *Scand J Gastroenterol* 2001; 36(3): 291-6.
  30. John B, Fagerhol MK, Lyberg T, Prydz H, Brandtzaeg P, Naess-Andresen CF, Dale I. Functional and clinical aspects of the myelomonocyte protein calprotectin. *Molecular Pathology* 1997; 50(3):113-23.
  31. Josefsson S, Bunn SK, Domellåf M. Faecal calprotectin in very low birth weight infants. *J Pediatr Gastroenterol Nutr.* 2007 Apr;44(4):407-13.
  32. Kristinsson J, Armbruster CH, Ugstad M, Kriwanek S, Nygaard K, Ton H, Fuglerud P. Faecal excretion of calprotectin in colorectal cancer; relationship to tumor characteristics. *Scand J Gastroenterol* 2001; 36(2):202-7.
  33. Kristinsson J, Roseth A, Fagerhol MK, Aadland E, Schjonsby H, Bormer OP, Raknerud N, Nygaard K. Faecal calprotectin concentration in patients with colorectal carcinoma. *Disease of the Colon & Rectum*1998; 41(3):316-21.
  34. Kronborg O, Ugstad M, Fuglerud P, John B, Hardcastle J, Scholefield JH, Vellacott K, Moshakis V, Reynolds JR. Faecal calprotectin levels in a high risk population for colorectal neoplasia. *Gut* 2000; 46(6):795-800.
  35. Kapel N, Campeotto F, Kalach N, Baldassarre M, Butel MJ, Dupont C. Faecal calprotectin in term and preterm neonates. *J Pediatr Gastroenterol Nutr.* 2010 Nov;51(5):542-7. Review.
  36. Kapel N, Roman C, Caldari D, Sieprath F, Canioni D, Khalifoun Y, Goulet O, Ruemmele FM. Faecal tumor necrosis factor-alpha and calprotectin as differential diagnostic markers for severe diarrhea of small infants. *J Pediatr Gastroenterol Nutr.* 2005 Oct;41(4):396-400.
  37. Kapel N, Barbot L, Gobert JG. New faecal markers: recent developments and perspectives. *Ann Pharm Fr.* 2004 Nov;62(6):371-5. Review. French.
  38. Kapel N. Faecal inflammatory markers in nutrition and digestive Disease in children. *Arch Pediatr.* 2004 May;11(5):403-5. Review. French. No abstract available.
  39. Laforgia N, Baldassarre ME, Pontrelli G, Indrio F, Altomare MA, Di Bitonto G, Mautone A. Calprotectin levels in meconium. *Acta Paediatr.* 2003 Apr;92(4):463-6.
  40. Limburg PJ, Ahlquist DA, Sandborn WJ, Mahoney DW, Devens ME, Harrington JJ, Zinsmeister AR. Faecal calprotectin levels predict colorectal inflammation among patients with chronic diarrhea referred for colonoscopy. *American Journal of Gastroenterology* 2000; 95(10):2831-7.
  41. Meucci G, D'Inca R, Maieron R, Orzes N, Vecchi M, Visentini D, Minoli G, Dal Pont E, Zilli M, Benedetti E, Virgilio T, Tonutti E. Diagnostic value of faecal calprotectin in unselected outpatients referred for colonoscopy: A multicenter prospective study. *Dig Liver Dis.* 2010 Mar;42(3):191-5. Epub 2009 Aug 19.
  42. Montalto M, Gallo A, Ferrulli A, Visca D, Campobasso E, Cardone S, D'Onofrio F, Santoro L, Covino M, Mirijello A, Leggio L, Gasbarrini G, Addolorato G. Faecal calprotectin concentrations in alcoholic patients: a longitudinal study. *Eur J Gastroenterol Hepatol.* 2011 Jan;23(1):76-80.
  43. Montalto M, Gallo A, Ianiro G, Santoro L, D'Onofrio F, Ricci R, Cammarota G, Covino M, Vastola M, Gasbarrini A, Gasbarrini G. Can chronic gastritis cause an increase in faecal calprotectin concentrations? *World J Gastroenterol.* 2010 Jul 21;16(27):3406-10
  44. Montalto M, Gallo A, Curigliano V, D'Onofrio F, Santoro L, Covino M, Dalvai S, Gasbarrini A, Gasbarrini G. Clinical trial: the effects of a probiotic mixture on non-steroidal anti-inflammatory drug enteropathy - a randomized, double-blind, cross-over, placebo-controlled study. *Aliment Pharmacol Ther.* 2010 Jul;32(2):209-14. Epub 2010 Apr 7.
  45. Montalto M, Santoro L, Dalvai S, Curigliano V, D'Onofrio F, Scarpellini E, Cammarota G, Panunzi S, Gallo A, Gasbarrini A, Gasbarrini G. Faecal calprotectin concentrations in patients with small intestinal bacterial overgrowth. *Dig Dis.* 2008;26(2):183-6. Epub 2008 Apr 21.
  46. Montalto M, Santoro L, Curigliano V, D'Onofrio F, Cammarota G, Panunzi S, Ricci R, Gallo A, Grieco A, Gasbarrini A, Gasbarrini G. Faecal calprotectin concentrations in untreated coeliac patients. *Scand J Gastroenterol.* 2007 Aug;42(8):957-61.
  47. Montalto M, Curigliano V, Santoro L, Armuzzi A, Cammarota G, Covino M, Mentella MC, Ancarani F, Manna R, Gasbarrini A, Gasbarrini G. Faecal calprotectin in first-degree relatives of patients with ulcerative colitis. *Am J Gastroenterol.* 2007 Jan;102(1):132-6. Epub 2006 Nov 13.
  48. Olafsdottir E, Aksnes L, Fluge G, Berstad A. Faecal calprotectin in infants with infantile colic, healthy infants, children with inflammatory bowel disease, children with recurrent abdominal pain and healthy children. *Acta Paediatr* 2002; 91: 45.
  - Roseth AG, Fagerhol MK, Aadland E, Schjonsby H. Assessment of the neutrophil dominating protein calprotectin in feces. A methodologic study. *Scandinavian Journal of Gastroenterology* 1992; 27(9):793-8.
  - Roseth AG, Aadland E, Jahnsen J, Raknerud N. Assessment of disease activity in ulcerative colitis by faecal calprotectin, a novel granulocyte marker protein. *Digestion*1997; 58(2):176-80.
  49. Roseth AG, Schmidt PN, Fagerhol MK. Correlation between faecal excretion of indium-111-labelled granulocytes and calprotectin, a granulocyte marker protein, in patients with inflammatory bowel disease. *Scandinavian Journal of Gastroenterology* 1999; 34(1):50-4.
  50. Roseth AG, Fagerhol MK, Aadland E, Schjonsby H. Assessment of the neutrophil dominating protein calprotectin in feces. A methodologic study. *Scand J Gastroenterol* 1992; 27(9):793-8.
  51. Roseth AG, Kristinsson J, Fagerhol MK, Schjonsby H, Aadland E, Nygaard K, Roald B. Faecal calprotectin: a novel test for the diagnosis of colorectal cancer? *Scand J Gastroenterol* 1993; 28(12):1073-6.
  52. Roseth AG, Aadland E, Grzyb K. Normalization of faecal calprotectin: a predictor of mucosal healing in patients with inflammatory bowel disease. *Scand J Gastroenterol.* 2004 Oct;39(10):1017-20.
  53. Scarpa M, D'Inca R, Basso D, Ruffolo C, Polese L, Bertin E, Luise A, Frego M, Plebani M, Sturmiolo GC, D'Amico DF, Angriman I. Faecal lactoferrin and calprotectin after ileocolonic resection for Crohn's disease. *Dis Colon Rectum.* 2007 Jun;50(6):861-9.
  54. Tibble JA, Bjarnason I. Department of Medicine, Guy's, King's, St Thomas's Medical School, Bessemer Road, London SE5 9PJ, UK. Non-invasive investigation of inflammatory bowel disease.
  55. Tibble JA, Bjarnason I. Department of Medicine, Guy's, King's, St. Thomas's Medical School, London, UK. Faecal calprotectin as an index of intestinal inflammation.
  56. Tibble JA, Sigthorsson G, Bridger S, Fagerhol MK, Bjarnason I. Surrogate markers of intestinal inflammation are predictive of relapse in patients with inflammatory bowel disease. [Journal Article] *Gastroenterology* 2000; 119(1):15-22.
  57. Tibble JA, Sigthorsson G, Foster R, Scott D, Fagerhol MK, Roseth A, Bjarnason I. High prevalence of NSAID enteropathy as shown by a simple faecal test. *Gut* 1999; 45(3):362-6.
  58. Tibble J, Teahon K, Thjodleifsson B, Roseth A, Sigthorsson G, Bridger S, Foster R, Sherwood R, Fagerhol M, Bjarnason I. A simple method for assessing intestinal inflammation in Crohn's disease. *Gut* 2000; 47(4):506-13.
  59. Tibble J, Sigthorsson G, Foster R, Fagerhol MK, Bjarnason I. Faecal calprotectin and faecal occult blood tests in the diagnosis of colorectal carcinoma and adenoma. *Gut* 2001; 49: 402.
  60. Tibble JA, Bjarnason I. Markers of intestinal inflammation and predictors of clinical relapse in patients with quiescent IBD. *Medscape Gastroenterol* 2001; 3 (2).
  61. Ton H, Brandsnes, Dale S, Holtlund J, Skuibina E, Schjonsby H, John B. Improved assay for faecal calprotectin. *Clinica Chimica Acta* 2000; 292(1-2):41-54.
  62. Kane, S., Sandborn, W., Rufo, P., Zhoulev, A., Boone, J., Lyster, D., Camilleri, M., and S. Hanauer. 2003. Faecal Lactoferrin is a Sensitive and Specific Marker in Identifying Intestinal Inflammation. *Am. J. Gastroenterol.* 98:1309-1314.
  63. Parsi, M., Shen, B., Achkar, J., Remzi, F., Goldblum, J., Boone, J., Lin, D., Connor, J., Fazio, V., and B. Lashner. 2004. Faecal Lactoferrin for the Diagnosis of Symptomatic Patients with Ileal Pouch-anal Anastomosis. *Gastroenterology* 126:1280-1286.
  64. Buderus, S., Boone, J., Lyster, D., and M. Lentze. 2003. Faecal Lactoferrin: New Parameter to Monitor Infliximab Therapy. *Dig. Dis. and Sciences.* 49:1036-1039.
  65. Langhorst, J., Elsenbruch, S., Mueller, T., Rueffer, A., Spahn, G., Michalsen, A., and G. Dobos. 2005. A Comparison Among Four Neutrophil-derived Proteins in Feces as an Indicator of Disease Activity in Ulcerative Colitis. *Inflamm. Bowel Dis.* 11(12):1085-1091.
  66. Kayazawa M, Saitoh O, Kojima K et al. Lactoferrin in whole gut lavage fluid as a marker for disease activity in inflammatory bowel disease: comparison with other neutrophil-derived proteins. *Am J Gastroenterol* 2002; 97: 360-369
  67. Tibble J, Teahon K, Thjodleifsson et al. A simple method for assessing intestinal inflammation in Crohn's disease. *Gut* 2000; 47: 506-513. van der Sluys Veer A, Brouwer J, Biemond I et al. Faecal lysozyme in assessment of disease activity in inflammatory bowel disease. *Dig Dis Sci* 1998; 43: 590-595
  68. Saitoh O, Sugi K, Matsue R et al. The forms and the levels of faecal PMN-elastase in patients with colorectal diseases. *Am J Gastroenterol* 1995; 90: 388-393
  69. Rachmilewitz D and International Study Group. Coated mesalazine versus sulfasalazine in the treatment of active ulcerative colitis: a randomized trial. *Br Med J* 1989; 298: 82-86
  70. WALKER C.W., "Faecal occult blood tests reduce colorectal cancer mortality.", *Am Fam Physician.* 2007 Jun 1;75(11):1652-3.
  71. CHIEN-HUA CHIANG, et al. «A comparative study of three faecal occult blood tests in upper gastrointestinal bleeding»; *Kaohsiung J. Med. Sci* May 2006, Vol 22, No 5: 223-228
  72. HIROFUMI MIYOSHI, et al. «Accuracy of Detection of Colorectal Neoplasia using an Immunohistochemical Occult Blood Test in Symptomatic Referred Patients: Comparison of Retrospective and Prospective Studies. *Internal Medicine* Sept. 2000 Vol. 39, No. 9: 701-706.
  73. Jemal A, Siegel R, Ward E, Hao Y, Xu J, Thun MJ. Cancer statistics, 2009. *CA Cancer J Clin* 2009; 59: 225-249
  74. Wei YS, Lu JC, Wang L, Lan P, Zhao HJ, Pan ZZ, Huang J, Wang JP. Risk factors for sporadic colorectal cancer in southern Chinese. *World J Gastroenterol* 2009; 15: 2526-2530

75. Lei T, Mao WM, Yang HJ, Chen XZ, Lei TH, Wang X, Ying Q, Chen WQ, Zhang SW. [Study on cancer incidence through the Cancer Registry Program in 11 Cities and Counties, China.]. *Zhonghua Liuxingbingxue Zazhi* 2009; 30: 1165-1170
76. Hewitson P, Glasziou P, Watson E, Towler B, Irwig L. Cochrane systematic review of colorectal cancer screening using the fecal occult blood test (hemoccult): an update. *Am J Gastroenterol* 2008; 103: 1541-1549
77. Logan RF. Review: faecal occult blood test screening reduces risk of colorectal cancer mortality. *Evid Based Med* 2009; 14: 15
78. Smith RA, Cokkinides V, Brooks D, Saslow D, Brawley OW. Cancer screening in the United States, 2010: a review of current American Cancer Society guidelines and issues in cancer screening. *CA Cancer J Clin* 2010; 60: 99-119
79. Mandel JS, Church TR, Bond JH, Ederer F, Geisser MS, Mongin SJ, Snover DC, Schuman LM. The effect of fecal occult-blood screening on the incidence of colorectal cancer. *N Engl J Med* 2000; 343: 1603-1607
80. Labianca R, Beretta GD, Kildani B, Milesi L, Merlin F, Mosconi S, Pessi MA, Prochilo T, Quadri A, Gatta G, de Braud F, Wils J. Colon cancer. *Crit Rev Oncol Hematol* 2010; 74: 106-133
81. Levin B, Brooks D, Smith RA, Stone A. Emerging technologies in screening for colorectal cancer: CT colonography, immunochemical fecal occult blood tests, and stool screening using molecular markers. *CA Cancer J Clin* 2003; 53: 44-55
82. Collins JF, Lieberman DA, Durbin TE, Weiss DG. Accuracy of screening for fecal occult blood on a single stool sample obtained by digital rectal examination: a comparison with recommended sampling practice. *Ann Intern Med* 2005; 142: 81-85
83. Allison JE, Sakoda LC, Levin TR, Tucker JP, Tekawa IS, Cuff T, Pauly MP, Shlager L, Palitz AM, Zhao WK, Schwartz JS, Ransohoff DF, Selby JV. Screening for colorectal neoplasms with new fecal occult blood tests: update on performance characteristics. *J Natl Cancer Inst* 2007; 99: 1462-1470
84. Oort FA, Terhaar Sive Droste JS, Van Der Hulst RW, Van Heukelem HA, Loffeld RJ, Wesdorp IC, Van Wanrooij RL, De Baaij L, Mutsaers ER, van der Reijt S, Coupe VM, Berkhof J, Bouman AA, Meijer GA, Mulder CJ. Colonoscopy-controlled intra-individual comparisons to screen relevant neoplasia: faecal immunochemical test vs. guaiac-based faecal occult blood test. *Aliment Pharmacol Ther* 2010; 31: 432-439
85. Parra-Blanco A, Gimeno-García AZ, Quintero E, Nicolás D, Moreno SG, Jiménez A, Hernández-Guerra M, CarrilloPalau M, Eishi Y, López-Bastida J. Diagnostic accuracy of immunochemical versus guaiac faecal occult blood tests for colorectal cancer screening. *J Gastroenterol* 2010; 45: 703-712
86. Young GP, Cole S. New stool screening tests for colorectal cancer. *Digestion* 2007; 76: 26-33
87. Kronborg O, Regula J. Population screening for colorectal cancer: advantages and drawbacks. *Dig Dis* 2007; 25: 270-273
88. Uchida K, Matsuse R, Miyachi N, Okuda S, Tomita S, Miyoshi H, Hirata I, Tsunoto S, Ohshiba S. Immunochemical detection of human blood in feces. *Clin Chim Acta* 1990; 189: 267-274
89. Burton RM, Landreth KS, Barrows GH, Jarrett DD, Songster CL. Appearance, properties, and origin of altered human hemoglobin in feces. *Lab Invest* 1976; 35: 111-115
90. Lengauer C, Kinzler KW, Vogelstein B. Genetic instability in colorectal cancers. *Nature* 1997; 386: 623-627
91. Tagore KS, Lawson MJ, Yucaitis JA, Gage R, Orr T, Shuber AP, Ross ME. Sensitivity and specificity of a stool DNA multitarget assay panel for the detection of advanced colorectal neoplasia. *Clin Colorectal Cancer* 2003; 3: 47-53
92. Woolf SH. A smarter strategy? Reflections on fecal DNA screening for colorectal cancer. *N Engl J Med* 2004; 351: 2755-2758
93. Hakama M, Coleman MP, Alexe DM, Auvinen A. Cancer screening: evidence and practice in Europe 2008. *Eur J Cancer* 2008; 44: 1404-1413
94. Hoff G, Grotmol T, Thiis-Evensen E, Bretthauer M, Gondal G, Vatn MH. Testing for faecal calprotectin (PhiCal) in the Norwegian Colorectal Cancer Prevention trial on flexible sigmoidoscopy screening: comparison with an immunochemical test for occult blood (FlexSure OBT). *Gut* 2004; 53: 1329-1333
95. McDonald S, Lyall P, Israel L, Coates R, Frizelle F. Why COMMENTS Chen JG et al. Colorectal cancer screening[WJG]www.wjgnet.com 2688 June 7, 2012[Volume 18]Issue 21[barium enemas fail to identify colorectal cancers. *ANZ J Surg* 2001; 71: 631-633
96. Kahi CJ, Imperiale TF, Juliar BE, Rex DK. Effect of screening colonoscopy on colorectal cancer incidence and mortality. *Clin Gastroenterol Hepatol* 2009; 7: 770-775; quiz 711
97. Rabeneck L, Paszat LF, Hilsden RJ, Saskin R, Leddin D, Grunfeld E, Wai E, Goldwasser M, Sutradhar R, Stukel TA. Bleeding and perforation after outpatient colonoscopy and their risk factors in usual clinical practice. *Gastroenterology* 2008; 135: 1899-1906, 1906.e1
98. Smith RA, Cokkinides V, Brawley OW. Cancer screening in the United States, 2008: a review of current American Cancer Society guidelines and cancer screening issues. *CA Cancer J Clin* 2008; 58: 161-179
99. Pox CP, Schmiegel W. Role of CT colonography in colorectal cancer screening: risks and benefits. *Gut* 2010; 59: 692-700
100. Chiang CH, Jeng JE, Wang WM, Jheng BH, Hsu WT, Chen BH. A comparative study of three fecal occult blood tests in upper gastrointestinal bleeding. *Kaohsiung J Med Sci* 2006; 22: 223-228
101. Sugi K, Saitoh O, Hirata I, Katsu K. Fecal lactoferrin as a marker for disease activity in inflammatory bowel disease: comparison with other neutrophil-derived proteins. *Am J Gastroenterol* 1996; 91: 927-934
102. Ward DG, Suggett N, Cheng Y, Wei W, Johnson H, Billingham LJ, Ismail T, Wakelam MJ, Johnson PJ, Martin A. Identification of serum biomarkers for colon cancer by proteomic analysis. *Br J Cancer* 2006; 94: 1898-1905
103. Ahmed N, Oliva KT, Barker G, Hoffmann P, Reeve S, Smith IA, Quinn MA, Rice GE. Proteomic tracking of serum protein isoforms as screening biomarkers of ovarian cancer. *Proteomics* 2005; 5: 4625-4636
104. Saitoh O, Kojima K, Kayazawa M, Sugi K, Tanaka S, Nakagawa K, Teranishi T, Matsuse R, Uchida K, Morikawa H, Hirata I, Katsu K. Comparison of tests for fecal lactoferrin and fecal occult blood for colorectal diseases: a prospective pilot study. *Intern Med* 2000; 39: 778-782
105. Hirata I, Hoshimoto M, Saito O, Kayazawa M, Nishikawa T, Murano M, Toshina K, Wang FY, Matsuse R. Usefulness of fecal lactoferrin and hemoglobin in diagnosis of colorectal diseases. *World J Gastroenterol* 2007; 13: 1569-1574
106. Sheng JQ, Li SR, Wu ZT, Xia CH, Wu X, Chen J, Rao J. Transferrin dipstick as a potential novel test for colon cancer screening: a comparative study with immuno fecal occult blood test. *Cancer Epidemiol Biomarkers Prev* 2009; 18: 2182-2185
107. Lönnnerdal B, Iyer S. Lactoferrin: molecular structure and biological function. *Annu Rev Nutr* 1995; 15: 93-110
108. Chew MH, Suzanah N, Ho KS, Lim JF, Ooi BS, Tang CL, Eu KW. Colorectal cancer mass screening event utilising quantitative faecal occult blood test. *Singapore Med J* 2009; 50: 348-353
109. Mandel JS, Church TR, Ederer F, Bond JH. Colorectal cancer mortality: effectiveness of biennial screening for fecal occult blood. *J Natl Cancer Inst* 1999; 91: 434-437
110. Allerberger F, et al. 2003. Hemolytic-uremic syndrome associated with enterohemorrhagic *Escherichia coli* O26:H infection and consumption of unpasteurized cow's milk. *Int. J. Infect. Dis.* 7:42-45. [PubMed]
111. Bastian S. N., Carle I., Grimont F. 1998. Comparison of 14 PCR systems for the detection and subtyping of stx genes in Shiga-toxin-producing *Escherichia coli*. *Res. Microbiol.* 149:457-472.[PubMed]
112. Bavaro M. F. 2009. *Escherichia coli* O157: what every internist and gastroenterologist should know. *Curr. Gastroenterol. Rep.* 11:301-306. [PubMed]
113. Belanger S. D., Boissinot M., Menard C., Picard F. J., Bergeron M. G. 2002. Rapid detection of Shiga toxin-producing bacteria in feces by multiplex PCR with molecular beacons on the Smart cyclor. *J. Clin. Microbiol.* 40:1436-1440. [PMC free article] [PubMed]
114. Bischoff C., Luthy J., Altwegg M., Baggi F. 2005. Rapid detection of diarrheagenic *E. coli* by real-time PCR. *J. Microbiol. Methods* 61:335-341. [PubMed]
115. 6. Chui L., et al. 2010. Comparison of Shiga toxin-producing *Escherichia coli* detection methods using clinical stool samples. *J. Mol. Diagn.* 12:469-475. [PMC free article] [PubMed]
116. 7. Cornu G., et al. 1999. Hemolytic uremic syndrome in Belgium: incidence and association with verocytotoxin-producing *Escherichia coli* infection. *Clin. Microbiol. Infect.* 5:16-22. [PubMed]
117. 119. Eklund M., Leino K., Siitonen A. 2002. Clinical *Escherichia coli* strains carrying stx genes: stx variants and stx-positive virulence profiles. *J. Clin. Microbiol.* 40:4585-4593. [PMC free article] [PubMed]
118. 120. Ethelberg S., et al. 2009. Outbreak of non-O157 Shiga toxin-producing *Escherichia coli* infection from consumption of beef sausage. *Clin. Infect. Dis.* 48:e78-81. [PubMed]
119. 121. Fey P. D., Wickert R. S., Rupp M. E., Safranek T. J., Hinrichs S. H. 2000. Prevalence of non-O157:H7 shiga toxin-producing *Escherichia coli* in diarrheal stool samples from Nebraska. *Emerg. Infect. Dis.* 6:530-533. [PMC free article] [PubMed]
120. 122. Gannon V. P., King R. K., Kim J. Y., Thomas E. J. 1992. Rapid and sensitive method for detection of Shiga-like toxin-producing *Escherichia coli* in ground beef using the polymerase chain reaction. *Appl. Environ. Microbiol.* 58:3809-3815. [PMC free article] [PubMed]
121. 123. Gilmour M. W., et al. 2009. Isolation and detection of Shiga toxin-producing *Escherichia coli* in clinical stool samples using conventional and molecular methods. *J. Med. Microbiol.* 58:905-911.[PubMed]
122. 124. Gould L. H., et al. 2009. Recommendations for diagnosis of shiga toxin-producing *Escherichia coli* infections by clinical laboratories. *MMWR Recomm. Rep.* 58:1-14.
123. 125. Griffin P. M. 1995. *Escherichia coli* O157:H7 and other enterohemorrhagic *Escherichia coli*, p. 739-761 In Blaser M. J., Smith P. D., Ravdin J. I., Greenberg H. B., Guerrant R. L., editors. (ed.), *Infections of the gastrointestinal tract*. Raven Press, New York, NY.
124. 126. Hanna S. S., Atkinson R. M., Jones T. F., Dunn J. R. 2010. Shiga toxin-producing *Escherichia coli* (STEC) in Tennessee: surveillance and current clinical laboratory practices. *Tenn. Med.* 103:39-42.[PubMed]
125. 127. Jelacic J. K., et al. 2003. Shiga toxin-producing *Escherichia coli* in Montana: bacterial genotypes and clinical profiles. *J. Infect. Dis.* 188:719-729. [PubMed]
126. 128. Johnson K. E., Thorpe C. M., Sears C. L. 2006. The emerging clinical importance of non-O157 Shiga toxin-producing *Escherichia coli*. *Clin. Infect. Dis.* 43:1587-1595. [PubMed]
127. 129. Karch H., Bielaszewska M. 2001. Sorbitol-fermenting Shiga toxin-producing *Escherichia coli* O157:H- strains: epidemiology, phenotypic and molecular characteristics, and microbiological diagnosis. *J. Clin. Microbiol.* 39:2043-2049. [PMC free article] [PubMed]
128. 130. Lockary V. M., Hudson R. F., Ball C. L. 2007. Shiga toxin-producing *Escherichia coli*, Idaho. *Emerg. Infect. Dis.* 13:1262-1264.
129. 131. Manning S. D., et al. 2007. Surveillance for Shiga toxin-producing *Escherichia coli*, Michigan, 2001-2005. *Emerg. Infect. Dis.* 13:318-321. [PMC free article] [PubMed]
130. 132. McCarthy T. A., et al. 2001. Hemolytic-uremic syndrome and *Escherichia coli* O121 at a lake in Connecticut, 1999. *Pediatrics* 108:E59. [PubMed]
131. 133. Misselwitz J., et al. 2003. Cluster of hemolytic-uremic syndrome caused by Shiga toxin-producing *Escherichia coli* O26:H11. *Pediatr. Infect. Dis. J.* 22:349-354. [PubMed]
132. 134. Nataro J. P., Kaper J. B. 1998. Diarrheagenic *Escherichia coli*. *Clin. Microbiol. Rev.* 11:142-201.[PMC free article] [PubMed]
133. 135. Nielsen E. M., Andersen M. T. 2003. Detection and characterization of verocytotoxin-producing *Escherichia coli* by automated 5' nuclease PCR assay. *J. Clin. Microbiol.* 41:2884-2893.[PMC free article] [PubMed]
134. 136. Paton A. W., Paton J. C. 2002. Direct detection and characterization of Shiga toxin-producing *Escherichia coli* by multiplex PCR for stx1, stx2, eae, ehxA, and saa. *J. Clin. Microbiol.* 40:271-274.[PMC free article] [PubMed]
135. 137. Schimmer B., et al. 2008. Outbreak of haemolytic uraemic syndrome in Norway caused by stx2-positive *Escherichia coli* O103:H25 traced to cured mutton sausages. *BMC Infect. Dis.* 8:41.[PMC free article] [PubMed]
136. 138. Spitzer M. D. 2002. Viral causes of diarrhea. *Pediatr. Rev.* 23:257-258. [PubMed]
137. 139. Thompson L. H., Giercke S., Beaudoin C., Woodward D., Wylie J. L. 2005. Enhanced surveillance of non-O157 verotoxin-producing *Escherichia coli* in human stool samples from Manitoba. *Can. J. Infect. Dis. Med. Microbiol.* 16:329-334. [PMC free article] [PubMed]
138. 140. Wang G., Clark C. G., Rodgers F. G. 2002. Detection in *Escherichia coli* of the genes encoding the major virulence factors, the genes defining the O157:H7 serotype, and components of the type 2 Shiga toxin family by multiplex PCR. *J. Clin. Microbiol.* 40:3613-3619. [PMC free article] [PubMed]
139. 141. Werber D., et al. 2002. A multistate outbreak of Shiga toxin-producing *Escherichia coli* O26:H11 infections in Germany, detected by molecular subtyping surveillance. *J. Infect. Dis.* 186:419-422.[PubMed]
140. 142. Wilhelmli I., Roman E., Sanchez-Fauquier A. 2003. Viruses causing gastroenteritis. *Clin. Microbiol. Infect.* 9:247-262. [PubMed]